

Thinky Application Datasheet 8



Carbon Nano Tube (1.0, 2.0 v %) THINKY Mixing Time - 6 minutes, 5Torr

Materials: (1) Carbon Nano Tube and (2) Epoxy polymer

Capacity: Total 60g

Material Ratio: Tube 1.0 to 2.0 v %

THINKY Mixer: Compare the Resistivity between Vacuum mixer ARV-310 and Vacuum-less ARE-250

THINKY mixing - 1v%													
	<p style="text-align: center;"><u>Resistivity</u></p>  <table border="1" style="margin-left: auto; margin-right: auto;"> <caption>Resistivity Data (k Ohm)</caption> <thead> <tr> <th>Mixing Method</th> <th>1 Vol %</th> <th>2 Vol %</th> </tr> </thead> <tbody> <tr> <td>Hand Mix</td> <td>~8,000</td> <td>~80</td> </tr> <tr> <td>ARE-250</td> <td>~100</td> <td>~20</td> </tr> <tr> <td>ARV-310</td> <td>~30</td> <td>~4</td> </tr> </tbody> </table>	Mixing Method	1 Vol %	2 Vol %	Hand Mix	~8,000	~80	ARE-250	~100	~20	ARV-310	~30	~4
Mixing Method	1 Vol %	2 Vol %											
Hand Mix	~8,000	~80											
ARE-250	~100	~20											
ARV-310	~30	~4											
<p>Nano Tube dispersion in viscous epoxy SEM Image(x20,000)</p>	<p>Vacuum deform has a lot of impact for resistivity improvement</p>												